

Nova-S RF Bipolaire

Application

1. Wrinkle removal, face lift and skin tighten, skin rejuvenation;
2. Redundant cellulite removal;
3. Body slimming and body shaping;
4. Fatigue relieving;
5. Channels and collaterals dredging;
6. Improve metabolism lymph detoxification.



Treatment Principle

Positioned tissues are heated to promote blood circulation and tighten the collagen in dermis. At the same time, the cooling measures are taken on the skin surface; the dermal layer of the epidermis is heated to maintain a normal temperature. Most of RF energy enter into the deeper dermal layer, so that deeper layers of skin to be up to 45 degree -60 degree, dermal layer of skin become thickened, and the wrinkles are shallow or disappear, so the skin become firmer and contour are lifted. The original collagen protein become thickened to increase the new collagen protein which are produced by irritate the skin.

Advantages

1. Most advanced technology.
2. New upgraded design and elegant appearance;
3. Suitable to all kinds of skin and people;
4. Comfortable, no pain, no scar in treatment;
5. High and adjustable power output, more convenient and effective;
6. Long working time: more than 10 hours without stopping per time;
7. Quick treatment, not affect people normal life;
8. 5.6 Inches color touch LCD screen of easy operation;
9. Unique external design, with high market competitiveness;
10. Multi functions, low cost and high profits.

Contraindication

1. Pregnant woman, abdomen (during menstrual period)
2. Audiphones, cardiac pacemaker, etc metal medical apparatus and steel nail, steel plate implant in the human body
3. Skin disease, the inflammation of the skin
4. Malignant tumors
5. Upper eyelid
6. Don't stay on the skin.

Treatment Item

1. Facial wrinkle removal (forehead wrinkle removal, remove wrinkle around eye, Horizontal wrinkle removal, removal wrinkle around mouse)
2. Neck wrinkle removal
3. Breast lifting
4. Body shaping(waist slimming, buttock lifting, leg slimming etc)

Technical Parameter

Bipolar RF energy	5%-100%(20 levels are available)
RF output	200W
Diameter of treatment head	25mm 45mm
Screen	5.6 Inches color touch screen
Work patterns	CW/PW
Net weight	10KG
Dimension	220mm×220mm×380mm(L×W×H)
Electrical requirement	AC220V±10% 50HZ AC110V±10% 60HZ



Interface

We designed the good operating system(5.6 inches Large color touch LCD screen), which makes more convenient operation and more exact treatment.



Professional treatment head

Small RF Head

wrinkle removal, face lift, skin tighten, improve blood circulation,etc

Φ25Head

face and neck lift,cellulite removal, improve metabolism

Φ45 Head

cellulite removal, body slimming, improve metabolismlymph detoxification



Treatment head



Big size treatment head: Apply to body.



Medium size treatment head: Apply to face, neck and hand.



Small size treatment head: Apply to rim of the eye, mouth and nose.

Treatment preparation

Conducting medium

Cooling Gel

Cooling gel is the ideal conducting medium. It could conduct heat more uniform. And the protect the skin. The whole treatment is comfortable under the help of cooling gel.



Cooling Gel



Distributeur exclusif

Nova-S

www.nova-s.fr